

Handheld *Retinomax Series 3*



Hand Held Ref
Retinomax 3




Handheld Autorefract Keratometer
Retinomax K-plus 3



The Retinomax Series, the world's leading

New Much lighter and has a thinner grip



Weighing just under one kilogram, the main body with battery is over 17% lighter in the case of the K-plus 3 and 11% lighter in the case of the Retinomax 3 than conventional Retinomax Series models.

The thinner grip makes it easier to hold



□ Retinomax 3
■ Retinomax 2



Retinomax 3

- Main body (with battery): 969 g
Down 11.1% from the conventional model (1,080 g)
- Retinomax K-plus 3
Main body (with battery): 999 g
Down 17.5% from the previous model (1,200 g)
- The girth of the grip (upper part): 145 mm
Down 11.5% from the previous model (175 mm)



The center of gravity is in the grip.

The Retinomax Series 3's center of gravity is in the grip, making the Retinomax Series 3 feel much lighter than it actually is.

New Auto pupil measurement, display and printout

The unit measures pupil size automatically, displays the X (horizontal) size on the monitor, and prints out the XY (horizontal and vertical) sizes separately. These results can be used for reference for checking accommodation, mydriasis, ADIE and Horner-associated syndrome, etc. The user can choose to print out either all the data or representative values only.



Pupil size displayed on the monitor

Representative values only	Whole values
Name: 107. 8. 7 3:48AM No. 018 VD: 13.75 -REF- [R] SPH CYL AX -2.25 -1.25 173 -2.25 -1.25 174 -2.25 -1.25 174 -2.25 -1.25 174 * -3.25 -2.25 174 10 x 8.1 y 8.1 [L] SPH CYL AX -2.25 -1.25 18 -2.25 -1.25 12 -2.25 -1.25 12 -2.25 -1.25 14 -2.25 -1.25 15 * -2.50 -1.25 14 8 x 5.7 y 5.4	Name: 107. 8. 7 3:52AM No. 028 VD: 13.75 -REF- [R] SPH CYL AX -2.75 -2.00 174 -2.75 -2.00 173 -2.75 -2.00 174 -2.75 -2.00 173 * -2.75 -2.00 173 10 x 8.5 y 8.5 [L] SPH CYL AX -2.50 -1.50 13 -2.25 -1.25 14 -2.25 -1.25 14 -2.00 -1.25 12 -2.00 -1.25 13 * -2.25 -1.25 13 8 x 5.7 y 5.5

Some data does not show pupil size due to eyelash obstruction or unstable instrument position.


New Auto Quick measurement

If no measurement is made within six seconds, the unit automatically enters Quick mode. Measurement is then taken in 0.07 second. In Auto Quick mode, "AQ" appears on printouts.


New Parallel sensor detects inclination of main body and displays value

Unlike conventional methods, in which operators guess levelness, levelness can be detected digitally with the use of the unit's parallel sensor. Levelness is displayed on the monitor in every 2° pitch up to ±12°, and every 3° pitch up to ±45°. (This is only displayed in the standing position. No levelness display appears when the patient is lying down.)

K-plus 3



Retinomax 3




Handheld Ref and Refract Keratometer, now boasts




New Alignment indicator display

Indicators showing alignment directions are displayed on the monitor to facilitate alignment. The Mire ring makes focusing much easier.



New Two memory functions

- 1 Data from the last measurement is automatically saved after the Retinomax 3 is turned off. The last patient data can be recalled by pressing the PRINT key. The data can be transferred wirelessly to a printer.
- 2 Press the MEMORY key for more than one second to enter Memory Set. Here you can save the data of up to 50 patients (100 eyes) in 10, 20, 30- and 40-patient steps. Data can also be printed out with Memory Set. For measurements, first press the MEMORY key to enter Memory mode. After each measurement, press PRINT to save individual patient data. This avoids the possibility of single-eye readings being mixed up. In normal readings, numbers are assigned in series. However, in Memory mode, saved data is numbered from "1."



New Achieving long and continuous use

With the adoption of a lithium ion battery, continuous operational time has been extended to approximately 80 minutes. The station has an AC adapter, allowing it to be connected to the main body with a DC cord (optional) to enable extended continuous operation. This facilitates screening, as well as use in an operating room.

DC cord: 2.9 m (length), 4 mm (diameter)



Optional DC cord connects main body and station

New Fixation intensity is automatically lowered for pupils less than 3 mm in size

The minimum pupil measurement size in Auto Quick mode is 2.3 mm. If a pupil is less than 3.0 mm, the fixation target intensity is automatically reduced by approximately 45% to avoid pupil contraction. The intensity can also be reduced with the fixation intensity key.

New Extended diopter adjustment range

The diopter adjustment range has been extended to $\pm 8D$.




Up to 50 sets of patient data can be saved or printed out. To delete saved data, select the DELETE ALL on the Memory Set screen.

To print out multiple sets of data, use the optional printer cable (main body to printer) to transfer data to the printer.



... greater mobility, stability and accuracy.

New Accuracy and stability have been improved



New Quick startup and fast printout

The main body takes only four seconds to start after power is turned on (some 60% faster than conventional models). Printout time has been reduced by 45%.

New Three fixation targets are available.

Specify when ordering.



```

'07. 8. 7 3:53AM
Name:
No. 027 VD: 13.75
-REF-
[R] SPH CYL AX
  -2.50 -2.00 176
  -2.50 -2.00 173
  -2.50 -2.00 174
  -2.50 -2.00 176
  -2.50 -2.25 178
* -2.50 -2.00 176 10
  x 5.7 y 6.7
[L] SPH CYL AX
  -2.00 -1.25 12
  -2.00 -1.50 8
  -2.00 -1.25 15
  -2.00 -1.25 9
  -2.00 -1.50 11
* -2.00 -1.25 11 10
  x 5.5 y 5.3
-KER-
[R] R1 R2 AX1 AX2
* 7.94 7.46 179 89
mm D deg
R1 7.94 42.50 179
R2 7.46 45.25 89
AV 7.70 43.87
CYL -2.75 179
[L] R1 R2 AX1 AX2
* 7.83 7.49 1 91
mm D deg
R1 7.83 43.12 1
R2 7.49 45.00 91
AV 7.66 44.00
CYL -1.88 1
    
```

```

'07. 8. 7 3:53AM
Name:
No. 027 VD: 12.0
[R] SPH CYL AX
  -2.50 -2.00 176
    
```

```

[L] SPH CYL AX
  -2.00 -1.25 11
    
```

New Retinomax 3 also has retro illumination mode

The Retinomax 3 includes the retro illumination mode for observation of the inside of the pupil. Accordingly, the only major difference between the Retinomax 3 and the Retinomax K-plus 3 is the adoption of kerato functions (including peripheral kerato readings).


In retro illumination mode



Retinomax 3



Startup screen



Measurement screen



Measurement completed screen

Retinomax K-plus 3



Startup screen



Measurement screen



Measurement completed screen


Pupil size: One
 Ref data: All
 Kerato data: Representative values

Eye diagram: ON

Order can be changed to either Right Ref & Kerato or Left Ref & Kerato

Stylish, streamlined designs and bright colors


Hand Held Ref
Retinomax 3



Operation panel (main body)

- 1 Fixation intensity key
- 2 MEMORY key
- 3 Axis rotation key
- 4 POWER key
- 5 QUICK/Melody key
- 6 RETRO mode/Pupil size key
- 7 R/L(right/left) select key
- 8 PRINT key
- 9 START switch

Handheld Autorefract Keratometer
Retinomax K-plus 3



Operation panel (main body)

- 1 Fixation intensity key
- 2 MEMORY key
- 3 PERI key (cornea peripheral measurement)
- 4 Axis rotation key
- 5 QUICK/Melody key
- 6 POWER key
- 7 MODE key
- 8 RETRO mode/Pupil size key
- 9 R/L(right/left) select key
- 10 PRINT key
- 11 START switch



Outstanding functionality in keeping with past Retinomax Series models

- The angle of the viewfinder can be adjusted within a range of 0° to 135°, making measurement easy regardless of patient's position or posture.



- A cylinder axis correction function with a 45° angle facilitates the use of the Retinomax on patients lying down.



- The 50-mm working distance is the same as that of the table-top model.
- The melody function provides a more relaxed atmosphere for children.



- Wireless data transfer from the main body to a printer is possible within a range of 8.7 meters directly and 6 meters at an angle of 30°.



Optional accessories



Hard carry case



Soft carry case



Spare battery



PV/PC interface cable




Printer cable
(main body + printer)



DC cord

Dimensions



Unit: mm

Specifications

Hand Held Ref Retinomax 3		
Refractometry		
Measurement range	Spherical (S + C)	-18.00D to +23.00D (in 0.25D increment)
	Cylinder	0 to +12D or 0 to -12D (in 0.25D increment)
	Cylinder axis	1 to 180° (in 1° increment)
Minimum pupil size	ø2.3 mm (Auto Quick)	
Vertex distance	0, 12 mm adjustable	
Measurement time	0.14 sec. per indication (continuous)	
Pupil size measurement time	+0.27 sec./1 time	
Measurement mode	Auto/Continuous/Auto Quick/Quick	
IDL wearing eye	Automatic support	
Fixation target	Picture target	
Battery life	Approx. 80 minutes (fully charged)	
Eye-piece	4.05 x 3.02 mm LCD	
Working distance	50 mm	
Reaching distance	278 mm	
External output	Infrared (Printer, Remote Vision)	
	RS-232C (Main body—PC, Printer—PC, Main body—Monitor)	
Dimensions (main body only)	170(W) x 230(H) x 240(D) mm	
Weight	Approx. 969 g (with battery) (Battery: Approx. 100 g)	

Station	
Dimensions	185(W) x 102(H) x 263(D) mm
Weight	Approx. 1.36 kg
Power consumption	100 VA

Printer	
Interface	RS-232C
Dimensions	93(W) x 77(H) x 263(D) mm
Weight	Approx. 595 g (without battery)

Handheld Autorefract Keratometer Retinomax K-plus 3		
Refractometry		
Measurement range	Spherical (S + C)	-18.00D to +23.00D (in 0.25D increment)
	Cylinder	0 to +12D or 0 to -12D (in 0.25D increment)
	Cylinder axis	1 to 180° (in 1° increment)
Minimum pupil size	ø2.3 mm (Auto Quick)	
Vertex distance	0, 12 mm adjustable	
Measurement of radius of curvature		
Measurement range	Radius of curvature	5.00 to 11.00 mm (in 0.01 mm increment)
	Corneal astigmatism	In 0.25D increment
	Cylinder axis	1 to 180° (in 1° increment)
	Center	ø3.2 mm (R8 mm)
	Peripheral	ø6.8 mm (R8 mm) (25° vertical/horizontal)
Measurement time	Kerato and Refract (continuous)	0.34 sec. per indication (continuous)
	Kerato alone	0.2 sec. per indication (continuous)
	Refract alone	0.14 sec. per indication (continuous), 0.07 sec. (Quick)
Pupil size measurement time	+0.27 sec./1 time	
Measurement mode	Auto/Continuous/Auto Quick/Quick	
IDL wearing eye	Automatic support	
Fixation target	Picture target	
Battery life	Approx. 80 minutes (fully charged, at 25°C)	
Eye-piece diopter adjustable range	±8D	
Working distance	50 mm	
Reaching distance	278 mm	
External output	Infrared (Printer, Remote Vision)	
	RS-232C (Main body—PC, Printer—PC, Main body—Monitor)	
Dimensions (main body only)	170(W) x 230(H) x 240(D) mm	
Weight	Approx. 999 g (with battery) (Battery alone: Approx. 100 g)	

Retinomax Series 3 meets CE (EMC, GM), UL and FDA standards

	WARNING: To ensure correct usage, read all manuals carefully before using equipment
--	--

Specifications and equipment are subject to change without any notice or obligation on the part of the manufacturer. © 2007/2008 RIGHT MFG. CO., LTD.
The information in this brochure is correct as of June 2008.

RIGHT MFG. CO., LTD.

Ophthalmic Sales

1-47-3, Maeno-cho, Itabashi-ku, Tokyo 174-8633, Japan
Tel: +81-3-3960-2275 Fax: +81-3-3960-2285
e-mail: eigyousitsu@rightmfg.co.jp

TOHOKU RIGHT MFG. CO., LTD.

Ophthalmic Service

45-1, Aza-yashikimae, Nakamura Osato-cho, Kurokawa-gun,
Miyagi 981-3521, Japan
Tel: +81-22-359-3113 Fax: +81-22-359-3413

Printed in Japan (0806-02)TI