

Sirius+ product improvements compared to the previous model include:

- USB 3.0 PC connectivity; improving performance and speeds up data the processing
- Double the resolution and image sampling providing outstanding high resolution Scheimpflug images;
- Allows up to 100 tomographic measurement sections per single acquisition
- Reduces the effect of hyper-scattering from sclera or from opaque corneal tissues
- Allows corneal imaging in two different fields of view with white light or fluoresceine
- A new central color camera, allows for a full evaluation of Dry Eye Disease. Providing:
 - Evaluation of symptoms
 - Non Invasive Tear film Break-up time (NiBut)
 - Meibomian glands analysis
 - Tear meniscus height measurement
 - Evaluation of conjunctival and fluorescein corneal staining, conjunctival and limbar hyperemia
 - Tearscope-like lipid layer evaluation

